



Granada

Model S743 | Class 3 Scooter

Granada

Model S743 | Class 3 Scooter



LCD dashboard.



Rear brake lights.



Easy access battery cover.



Large front basket with handle.



Front LED light.



Comfortable height adjustable, reclining seat.



York Road
Bridgend Industrial Estate
Bridgend
CF31 3TB

Telephone: 01656 674488
Facsimile: 01656 674499
Email: sales@romamedical.co.uk

www.romamedical.co.uk

*Lower limb amputees should refer to the user manual before use.
** Passenger weight, ground conditions and battery condition/size can affect distance of travel.
*** Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.

Special Features

- Attractive removable modesty panel.
- Full LED lights set.
- LCD dashboard for greater clarity.
- Dual mode wig-wag lever.
- Automatic speed reduction. (On full lock turns in 8mph setting for safety)
- Flip up, width and height adjustable armrests.
- Infinite angle tiller adjustment.
- Rear disc brake system.
- Non-marking pneumatic tyres.
- Comfortable height adjustable, reclining seat.
- Large front basket with carry handle.
- Automatic brake lights.
- Full adjustable suspension.

Specifications

Model No.*:	S743
Class (Invalid Carriage):	3
Overall Length:	1300 mm / 51"
Overall Width:	660 mm / 26"
Overall Height:	640 mm / 25"
Seat Width:	1240 mm / 49"
Safe Working Slope:	8° / 14%
Maximum User Weight:	160 kg / 25 st
Battery Size:	12v x2 / 50 Ah
Range**:	40 Km / 25 Miles
Overall Weight:	110 Kg / 242 lbs
Maximum Speed:	6 & 12 kmph / 4 & 8 mph
Wheel Size:	300 mm / 12"
Turning Radius:	1265 mm / 50"
Battery Charger:	6 Amp Standard
Ground Clearance***:	120 mm / 12"
Colours:	Red

Dealer Stamp

