

SHOPPRIDER
The World's Leading Scooter Manufacturer



Shoprider Grande

Class 3 Scooter

ROMA
MEDICAL

ROMA

SHOPRIDER® GRANDE



Special Features

- Full LED lights set
- Full suspension
- 50 amp power battery
- High and low speed settings
- Fully adjustable captain seat
- Black non-marking wide profile pneumatic tyres
- Free wheel release handle
- Delta handle
- USB charging port
- Flip up adjustable armrests
- Tiller angle adjustment lever
- Rear indicators & hazard lights

Specifications

| | | | |
|----------------------------------|-----------------------|---------------------------|-----------------------------|
| Model No.* | TE-9AS | Safe Working Slope | 8° / 14% |
| Class (Invalid Carriage) | 3 | Maximum User Weight | 159 kg / 25 st |
| Overall Length | 1400 mm / 55" | Battery Size 12v x 2 | 36Ah / 50Ah (Option) |
| Overall Width | 600 mm / 24" | Range** | 28 km / 18 miles |
| Overall Height (folded, no seat) | 1360 mm / 54" | Weight | 99 kg / 218 lbs |
| Overall Height (top of mirror) | 1170 mm / 46" | Maximum Speed | 12 kmph / 7.5 mph |
| Seat Width | 460mm / 18" | Wheel / Tyre Size (Front) | 330 mm / 13" |
| Seat Depth | 430 mm / 17" | (Back) | 330 mm / 13" |
| Seat Height (to ground) | 660 mm / 26" | Turning Radius | 1680 mm / 66" |
| Backrest Height (from seat) | 730 mm / 29" | Battery Charger | 5 amp standard |
| Seat Height (to floorpan) | 470 mm / 18.5" | Ground Clearance*** | 10mm / 4" |

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.
*** Measurement to underside of chassis



ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB
Telephone : 01656 674488 Facsimile : 01656 674499
Email : sales@romamedical.co.uk

www.romamedical.co.uk

Dealer Stamp