

ROMA
MEDICAL



Roma Marbella

Class 2 Powerchair

ROMA
MEDICAL

ROMA MARBELLA



Adjustable angle seat and backrest



Easily folded for storage and transportation



Removable battery pack



Swing in and out detachable footrests

Special Features

- Adjustable angle seat and backrest
- Adjustable height and width arms
- Twin detachable battery packs
- Backrest folds down for transport and storage
- Heaviest part only 18 kgs (main chassis)
- Heaviest part only 18 kgs (main chassis)
- Rear suspension
- 2" standard cushion
- Lap belt
- Optional elevating legrests
- Controller can be fitted on both sides
- Crash tested

Specifications

Model No.*	P200	Maximum User Weight	136 kg / 21 st
Class	2	Battery Size	2 x 12v 50 amp
Overall Length	1066 mm / 42"	Charger	5 amp off board
Overall Width	610 mm / 24"	Range**	40 km / 25 miles
Overall Height	960 mm / 38"	Weight	80 kgs / 176 lbs
Seat Width	460 mm / 18"	Maximum Speed	6 kmph / 4 mph
Seat Depth	400 mm / 16"	Wheel / Tyre Size	360 mm / 14"
Seat Height	495 - 520 mm / 19" - 20"	Turning Radius	580 mm / 23"
Backrest Height (from seat)	430 mm / 17"	Ground Clearance	100 mm / 4"
Safe Working Slope	8° / 12%	Colour	Red

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.

ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB
Telephone : 01656 674488 Facsimile : 01656 674499
Email : sales@romamedical.co.uk

www.romamedical.co.uk

Dealer Stamp

