

**ROMA**

# *WHEELCHAIR POWER PACK*



**ROMA** *Global Independence*

# ROMA WHEELCHAIR POWER PACK



## Special Features

- Twin drive wheels for greater traction
- Dismantles easily for easy transportation and storage
- Easy to operate hand control lever
- Variable speed setting with reverse
- Security key switch
- Removeable battery pack for remote charging
- Lifts out of use to enable normal wheelchair function
- Suitable for users up to 136 kgs / 21 stone



Mounting Brackets



Dismantles for easy transportation and storage



Power pack controller



Power pack

## Specifications

Model No.	?????	Maximum Range	10 miles / 16 kms
Battery Capacity	12v / 20Ah	Fits Wheelchair Width	360-510 mm / 14-20"
Charger (off board)	1.5A	Total Weight	13.2 kgs / 29 lbs
Drive Wheels (puncture proof)	2 x 200 mm / 8"	Motor Weight	6.4 kgs / 14 lbs
Maximum Drive Speed	3.4 mph / 5.5 kmph	Battery Weight	6.8 kgs / 15 lbs

CE

# ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB  
 Telephone : 01656 674488 Facsimile : 01656 674499  
 Email : sales@romamedical.co.uk

[www.romamedical.co.uk](http://www.romamedical.co.uk)

Dealer Stamp