



# Roma Valencia

Class 2 Scooter

## Accessories



Oxygen bottle holder



Rear basket & attachment



Crutch & walking stick holder



Canopy



Rear storage box

RRP £2,499.00 Our Price £

# ROMA VALENCIA



Easily dismantles for storage and transportation



Full lights set including front indicators



Comfortable reclining seat with slide and swivel base



Split battery pack

## Special Features

- Comfortable seat - lockable four position swivel
- Height adjustable seat
- Width adjustable arms
- Headlight and rear reflectors
- Adjustable angle tiller
- Split battery pack
- Dismantles easily for storage and transportation
- Black non-marking pneumatic tyred wheels for extra comfort

## Specifications

Model No.*	<b>GK10</b>	Safe Working Slope	<b>8° / 14%</b>
Class	<b>2</b>	Maximum User Weight	<b>135 kg / 21 st</b>
Overall Length	<b>1170 mm / 46"</b>	Battery Size 12v x 2	<b>34/36Ah</b>
Overall Width	<b>585 mm / 23"</b>	Range**	<b>32 km / 20 miles</b>
Overall Height	<b>960 mm / 38"</b>	Weight (with Battery Pack, Seat, Basket)	<b>75 kgs / 135 lbs</b>
Seat Width	<b>460 mm / 18"</b>	Maximum Speed	<b>6 kmph / 4 mph</b>
Width Between Arms	<b>460 - 610 mm / 18" - 24"</b>	Wheel / Tyre Size	<b>250 mm / 10"</b>
Seat Height (to ground)	<b>580 - 655 mm / 23" - 26"</b>	Turning Radius	<b>1330 mm / 52"</b>
Armrest Height	<b>135 - 210 mm / 5" - 8"</b>	Battery Charger	<b>4 amp off board</b>
Backrest Height from Seat	<b>440 mm / 16"</b>	Ground Clearance***	<b>90 mm / 3.5"</b>
Seat to Footrest/Floorpan Height	<b>480mm - 555mm / 19" - 22"</b>	Colours	<b>Metallic Red / Blue</b>

\* Lower limb amputees should refer to user manual before use. \*\* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.  
\*\*\* Measurement to underside of chassis



# ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB  
Telephone : 01656 674488 Facsimile : 01656 674499  
Email : sales@romamedical.co.uk

[www.romamedical.co.uk](http://www.romamedical.co.uk)

Dealer Stamp