



# Sorrento

Model S742 | Class 2 Scooter

# Sorrento

Model S742 | Class 2 Scooter



Non marking puncture proof tyres.



Fixed Frame.



Attractive front and rear LED lights.



Front shopping basket with carry handle.



Infinite angle tiller adjustment.



Flip up and width adjustable arm rests.

## Special Features

- Dual charge sockets for on and off board charging.
- Fixed frame.
- Non marking puncture proof tyres.
- Comfortable adjustable height seat.
- Flip up and width adjustable arm rests.
- Infinite angle tiller adjustment.
- Attractive front and rear LED lights.
- Front shopping basket with carry handle.

## Specifications

Model No.*:	S742
Class (Invalid Carriage):	2
Overall Length:	1200mm / 47"
Overall Width:	580mm / 23"
Overall Height:	1000mm / 39.5"
Seat Width:	460mm / 18"
Safe Working Slope:	8° / 14%
Maximum User Weight:	135kg / 300lbs
Battery Size:	2 x 34 / 36Ah
Range**:	Up to 32 km / 20 Miles
Overall Weight:	85kgs / 187lbs
Maximum Speed:	6 kmph / 4 mph
Wheel Size (Front):	230mm / 9"
Wheel Size (Back):	255mm / 10"
Turning Radius:	1430mm / 56"
Battery Charger:	5 Amp off board
Ground Clearance***:	90mm / 3.5"
Colours:	Red



York Road  
Bridgend Industrial Estate  
Bridgend  
CF31 3TB

Tel: 01656 674488  
Fax: 01656 674499  
Email: sales@romamedical.co.uk

**www.romamedical.co.uk**

 We are social. Follow us on Facebook  
/RomaMedicalAids

\*Lower limb amputees should refer to the user manual before use.  
\*\* Passenger weight, ground conditions and battery condition/size can affect distance of travel.  
\*\*\* Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.



Dealer Stamp