



Tulsa

Model S245 | Class 2 Scooter

Tulsa

Model S245 | Class 2 Scooter



Available in Orange & Red.



Front & rear LED headlight.



Easy to use controls.



Front & rear shock absorbers.



Wing mirrors.



Comfortable, height adjustable swivel seat.



York Road
Bridgend Industrial Estate
Bridgend
CF31 3TB

Tel: 01656 674488
Fax: 01656 674499
Email: sales@romamedical.co.uk

www.romamedical.co.uk

 We are social. Follow us on Facebook
[/RomaMedicalAids](https://www.facebook.com/RomaMedicalAids)

*Lower limb amputees should refer to the user manual before use.
** Passenger weight, ground conditions and battery condition/size can affect distance of travel.
*** Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.



Special Features

- Comfortable height & depth adjustable swivel seat.
- Front & rear LED lights.
- Front shopping basket.
- Wing mirrors.
- Shopping basket.
- Easy to use controls.
- Front & rear shock absorbers.
- Pneumatic tyres.
- Flip up arm rests.
- Padded and fold down backrest.
- Angle adjustable tiller.
- USB Charging port.

Specifications

Model No.*:	S245
Class (Invalid Carriage):	2
Overall Length:	1300mm / 51.2"
Overall Width:	610mm / 24"
Overall Height:	1170mm / 46"
Seat Width:	480mm / 18.8"
Safe Working Slope:	8°
Maximum User Weight:	136 kg / 21st
Battery Size:	38 Ah
Range**:	38 km / 24 miles
Overall Weight:	94 kg
Maximum Speed:	6 kmph (4 mph)
Wheel Size Front:	228.8mm / 9"
Wheel Size Rear:	254mm / 10"
Turning Radius:	1550mm / 61"
Battery Charger:	3 Ah
Ground Clearance***:	110mm / 4.3"
Colours:	Orange, Red

Dealer Stamp