

# Paediatric Self-Propelling Wheelchair

## Model No. 1451 / 1451 EL

- ▶ Vertical swing away, height adjustable footrests
- ▶ 56 cm / 22" pneumatic tyre, rear Mag wheels with steel hand rim
- ▶ 15 cm / 6" front solid tyre swivel castors
- ▶ Designed for children / young adults
- ▶ Elevating leg rests optional (AC103)
- ▶ Chrome plated steel frame
- ▶ Crash tested
- ▶ Progressive and parking braking system
- ▶ Nylon cord upholstered seat and back - flame retardant



## Dimensions

Depth	95 - 104cm / 37 - 41"
Width	56cm / 22"
Height	93cm / 36"
Seat height	46cm / 18"
Seat depth	33cm / 13"
Seat width	35.5cm / 8-11"
Seat to foot rest	30-39 cm / 12 - 15"
Elevating legrests	20 -28 cm / 8 - 11"
Weight of product	15kgs / 33lbs
Maximum user weight	83 kgs / 183 lbs

# Child's Rollator

## Model No. 2243

- ▶ Chrome plated steel
- ▶ Adjustable handles
- ▶ Solid tyre wheels
- ▶ Rear non-slip rubber feet

## Dimensions

Depth	47 cm / 19"
Maximum width	42 cm / 17"
Height	97-142 cm / 38-56"
Weight of product	6 kgs / 13 lbs
Maximum user weight	75 kgs / 165 lbs



## Paediatric Adjustable Crutch

### Model No. 2120

- ▶ Open arm band
- ▶ Lightweight anodised tubing
- ▶ Non-slip crutch tip
- ▶ Moulded handgrip and cuff

### Dimensions

Cuff to grip accordance with 18 cm / 7"  
EN ISO 11334-17

Grip to ground in accordance 56cm / 22"  
with EN ISO 11334-124



## Chrome Paediatric Adjustable Tripod

### Model No. 2474/AM

- ▶ Swan neck Chrome plated steel handle
- ▶ Plastic formed grip
- ▶ Chrome plated steel base
- ▶ Moulded polypropylene seat
- ▶ Non-slip, non-marking rubber tips
- ▶ Right or left hand use

### Dimensions

Height 48 - 58cm / 19 - 23"

Width 30cm / 12"

Depth 20cm / 8"



## Coloured Paediatric Adjustable Tripod

### Model No. 2516

- ▶ Swan neck epoxy coated steel handle
- ▶ Plastic formed grip
- ▶ Epoxy coated steel base
- ▶ Moulded polypropylene seat
- ▶ Non-slip, non-marking rubber tips
- ▶ Right or left hand use

### Dimensions

Height 48 - 58cm / 19 - 23"

Width 30cm / 12"

Depth 20cm / 8"



# Child's Rollator

## Model No. 2243

- ▶ Chrome plated steel
- ▶ Adjustable handles
- ▶ Solid tyre wheels
- ▶ Rear non-slip rubber feet

### Dimensions

Depth	47 cm / 18.5"
Maximum width	42 cm / 16.5"
Height	42-56cm / 16.5 - 22"
Weight of product	6 kgs / 13 lbs
Maximum user weight	75 kgs / 165 lbs

