

Commodes & Sanichairs

Romachair
Mobile Commodes
Stacking Commodes
Commode Accessories





Romachair

Anti-Tip Hospital Commode

Several years ago, we were approached by numerous healthcare professionals with a critical need: the development of a hospital commode/toilet chair designed to effectively mitigate the spread of infections such as MRSA and Clostridium Difficile.

This led to the conception of the Romachair, meticulously engineered to meet stringent criteria: functionality, durability, ease of cleaning, adaptability for bedside use, and compatibility with existing hospital ward toilets.

Distinguished from its counterparts, the entire Romachair range is proudly manufactured in the UK, adhering rigorously to the highest British and international standards. Backed by our commitment to quality, we offer a comprehensive 5-year limited warranty on the frame, ensuring hospitals receive unparalleled value.

Featuring removable commode seats, over seats, and backrests, cleaning is effortlessly facilitated without the need to use any specialised tools.

Our readily available replacement parts ensure the longevity of the commode, sustaining its serviceability for many years to come.

RMA Health is ISO 9001, ISO-14001 and 13485 accredited using only the best quality material and components in the manufacturing of the Romachair.



Designed and Manufactured
in Great Britain

Romachair

Romachair Hospital commodes are available in either an anti microbial plastic coated steel finish or in anti microbial chemically polished stainless steel.

All of our Marine Grade anti microbial Stainless Steel (SS) Romachairs are hygienically cleaned after manufacture by a process known as Electro Chemical Polishing and Passivation. This process is extensively used throughout the pharmaceutical and medical industry. It renders the surface of the stainless steel up to 30 times more corrosive resistant than untreated stainless steel.

Model shown
4740SS (Stainless Steel)



All of our Electropolished Stainless Steel (SS) Romachairs are hygienically cleaned after manufacturing by a process known as electro chemical polishing and passivation.

Standard Romachair 4740

The 4740 Series covers patients weighing up to 250kg (40 Stone). Medium Romachair with a 21" (53cm) seat width, detachable arms and heel support.* Available in an Anti-Microbial powder coating or a chemically polished stainless steel finish.


Bariatric Romachair 4750A

The 4750 Series covers patients weighing up to 318kg (50 Stone). Bariatric Romachair with a 25" (63.5cm) seat width, fixed arms and heel support. Available in an Anti-Microbial powder coating or a chemically polished stainless steel finish.

Standard



CV Mod-



In hospital environments, two main types of toilet configurations are commonly encountered: toilets with wall-mounted cisterns and toilets where the cistern is either concealed behind a panel or integrated into the wall.


RMA Health received inquiries from several hospitals with toilets featuring concealed cistern designs. They reported issues arising from their current commode chairs, as the aperture did not align properly with the toilet bowl. This misalignment posed significant challenges for nursing and cleaning staff.

In response, RMA endeavoured to design a commode chair specifically tailored to resolve these issues and ensure seamless alignment with the toilet, resulting in the development of the Romachair CV model.

Seat Options

We have a range of seat options to suit all needs including the gap front seat, the full aperture seat. All of our commodes come supplied with a grey overseat cover, full aperture seat and pulp pan rack as standard.

Optional Blue Aperture Seat & Backrest Available



Commode Model	Gap Front	Full Aperture	Overseat Cover	Back Rest
4740	PL0218/G	PL0217/G	PL0235	PL0215/G
4750	NA	PL0216/G	PL0236	PL0215/G

Pan Racks

RMA Health also has the most extensive range of pan racks available in the market and when combined with any of our chairs we are able to cover almost any eventuality. There is no standard specification for bed pans in the UK. Many hospitals use the disposable pulp version but there are several other types in daily use. We mention below the majority of those now in use. (Pulp pan rack is supplied unless requested otherwise.)



*All pan racks are free floating, eliminating screws and attachments and can be decontaminated in a bed pan washer.

It is crucial to accurately specify the type of pan, aperture seat, and seat cover needed for the commode chair when placing an order to ensure the correct manufacturing and supply of the product.



	Vector Order No.	S/Steel Order No.	Warwick/Sasco Order No.	Orange Order No.	Square Order No.	Pulp Order No.	Blue Henley Order No.	ETAC Order No.
4740	AC163/M	AC151/M	AC157/M	AC166/M	AC162/M	AC158/M	AC165/M	AC168/M
4750	AC163/B	AC151/B	AC157/B	AC166/B	AC162/B	AC158/M		

Please Note: Commodes are supplied with pulp pan racks as standard unless requested otherwise.



Medium, Detachable
21" (53cm) Seat



Detachable Arms



Intergrated Anti-Tip
Heel Support



Weight capacity of
up to 250kg (40st)



4740 - 4740/CV
Antimicrobial Plastic
Coated



4740/SS - 4740/CV/SS
Chemically Polished
Stainless Steel

Romachair 4740

Model shown 4740

NHSSC - FYB 524



Model shown 4740SS

NHSSC - FYB 920



- Sealed Seamless Frame
- No Hinges on Foot/Heel Plate
- Suitable for all Commode Pans
- Comes complete with over-seat cover.
- Optional gap or full aperture seat available.
- Clipless Back & Seat

Romachair 4740/CV Straight Back

Model shown 4740CV

NHSSC - GTB 1408



Straight Back for
Concealed Cistern

Model shown 4740SS
with optional blue seat and back

NHSSC - GTB 1415



- Sealed Seamless Frame
- No Hinges on Foot/Heel Plate
- Suitable for all Commode Pans
- Comes complete with over-seat cover.
- Optional gap or full aperture seat available.
- Clipless Back & Seat

Model Code	Product Weight	Overall Height	Overall Width	Overall Depth	Seat Pad Height	Seat Pad Depth	Width Between Arms	Clearance to Under Seat	Max. User Weight
4740	18kg (34lbs)	104cm (41")	61cm (24")	66cm (26")	57cm (22.5")	45cm (18")	52cm (20.5")	51cm (20")	250kg (40st)
4740/CV	18kg (34lbs)	102cm (40")	61cm (24")	73cm (29")	52cm (20.5")	45cm (18")	52cm (20.5")	46cm (18")	250kg (40st)

Romachair 4750A

Model shown 4740/SS

NHSSC - GTB 1419

Model shown 4750A

NHSSC - GTV 85337



Bariatric, Detachable
25" (53cm) Full



Fixed Arms



Integrated Anti-Tip
Heel Support

Romachair 4750/CV Straight Back

Straight Back for
Concealed Cistern -

Model shown 4750/CV
with optional blue seat & back

NHSSC - GTV 85333



Model shown 4750/CV/SS

NHSSC - GTV 85345



Weight capacity of
up to 318kg (50st)



4750A - 4750/CV
Antimicrobial Plastic
Coated



4750/SS - 4750/CV/SS
Chemically Polished
Stainless Steel

Model Code	Product Weight	Overall Height	Overall Width	Overall Depth	Seat Pad Height	Seat Pad Depth	Width Between Arms	Clearance to Under Seat	Max. User Weight
4750A	18kg (39lbs)	104cm (41")	74.5cm (29")	81cm (32")	58.5cm (23")	46cm (18")	67cm (26")	51cm (20")	318kg (50st)
4750/CV	18kg (39lbs)	102cm (40")	71cm (28")	90cm (35.5")	52cm (20.5")	46cm (18")	63cm (25")	46cm (18")	318kg (50st)

Mobile Commode

Height Adjustable Drop Arm

- Chrome plated steel
- 4 brake castors 12.5 cm / 5"
- Drop down padded armrests
- Round Bowl - Slides to remove
- Flame retardant upholstered over seat and back rest
- Height adjustable legs

Dimensions

Overall Height	103cm (40.5")
Overall Width	53cm (21")
Overall Depth	58cm (23")
Seat Height	53cm (21")
Seat Depth	48cm (19")
Width between arms	45cm (17.7")
Toilet Clearance	43cm (17")

3190A/4BC

Optional Model no. AC01

Perfection pan rack

Optional Model no. PL0068

Grey round bowl

Mobile Commode

Heavy Duty Chair

- Chrome plated steel
- Detachable armrests
- Aperture seat
- Round Bowl - Slides to remove
- Flame retardant upholstered over seat and backrest
- Push handle
- 4 brake castors 12.5 cms / 5"

Dimensions

Overall Height	98cm (38.5")
Overall Width	63.5cm (25")
Overall Depth	56cm (22")
Seat Height	47cm (18.5")
Seat Depth	46cm (18")
Width between arms	53cm (21")
Toilet Clearance	50cm (19")
Maximum user weight	190kg (30st)

3275/4BC

Optional Model no. 3175/FR

Hook on footrests

Optional Model no. PL0068

Grey round bowl



MODEL NO. 3190/FR
Push on footrests Suitable for Model no. 3190



MODEL NO. 3175/FR
Hook on footrests Suitable for Model Nos. 3175, 3275 and 3375

Mobile Commode Detachable Backrest

- Chrome plated steel
- Detachable arms
- Aperture Seat
- Round Bowl - Slides to remove
- Flame retardant upholstered over seat and backrest
- Push handles

Dimensions

Overall Height	99cm (39")
Overall Width	72cm (28")
Overall Depth	63.5cm (25")
Seat Height	47cm (18.5")
Seat Depth	48cm (19")
Width between arms	58.5cm (23")
Toilet Clearance	50cm (19")
Maximum user weight	318kg (50st)



3375/4BC

Optional Model no. 3175/FR
Hook on footrests

Optional Model no. PL0068
Grey round bowl

Mobile Commode Detachable Backrest

- 12.5cm swivel type brake castors
- Detachable back
- Push handle
- Flame retardant upholstered over seat and backrest
- Round bowl – removable
- Toilet clearance is 50cm

Dimensions

Overall Height	98cm (38.5")
Overall Width	52cm (20")
Overall Depth	58cm (23")
Seat Height	55cm (21.5")
Seat Depth	47cm (18.5")
Width between arms	45cm (17.7")
Toilet Clearance	50cm (19")
Maximum user weight	155kg (25st)



3175/4BC

Optional Model no. 3175/FR
Hook on footrests

Optional Model no. PL0068
Grey round bowl



MODEL NO. AC01
Perfection pan rack Suitable for Model Nos. 3175, 3275, 3375 and 3190



MODEL NO. PL0068
Round Bowl Suitable for Model Nos. 3175, 3275, 3375 and 3190

Chemical Commode Chair

- White plastic coated steel
- Height adjustable legs
- Moulded armrests
- Lift up toilet seat with lid
- 20 litre plastic pan (removable) for chemical toilet liquid

Dimensions

Seat height	47-61 cm / 19-24"
Seat width	38 cm / 15"
Seat depth	43 cm / 17"
Width	59 cm / 23"
Height	71-86 cm / 28-34"
Depth	63.5 cm / 25"
Width between arms	46 cm / 18"
Product Weight	7 kgs / 15 lbs
Maximum user weight	160 kgs / 352 lbs

3308



Elsan Toilet Stand

- For use with Elsan Bristol Chemical toilet
- White plastic coated steel
- Plastic armrests
- Height adjustable legs

Dimensions

Seat height	54-64 cm / 21-25"
Seat width	44 cm / 17"
Seat depth	47 cm / 19"
Width	54.5 cm / 21.5"
Height	89 cm / 35"
Depth	58 cm / 23"
Width between arms	43 cm / 17"
Maximum user weight	155 kgs / 341 lbs
Weight of toilet	13.1 kgs / 29 lbs

3325



Optional Model no. 3326

Elsan Bristol toilet only

Optional Model no. 3327

Elsan blue fluid 8x2 litre case



Note: Toilet & Stand are sold separately.

Stacking Commode Armchair

- Traditional armchair appearance
- Plastic coated steel
- Hygienic pot (removable) with lid
- Upholstered backrest and overseat - flame retardant

Dimensions

Seat height	51 cm / 20"
Seat height (with overseat)	54 cm / 21"
Seat width	40 cm / 16"
Seat depth	40 cm / 16"
Width	54.5 cm / 21.5"
Height	84 cm / 33"
Depth	53 cm / 21"
Width between arms	49.5 cm / 19.5"
Weight of product	6.9 kgs / 15 lbs
Maximum user weight	133 kgs / 292 lbs

Optional Model no. 3510A
Adjustable Height



Stacking Commode

- Plastic coated steel
- Hygienic pot (removable) with lid
- Upholstered backrest and overseat - flame retardant

Dimensions

Seat height	49.5 cm / 19.5"
Seat height (with overseat)	54.5 cm / 21.5"
Seat width	43 cm / 17"
Seat depth	43 cm / 17"
Width	56 cm / 22"
Height	89 cm / 35"
Depth	53 cm / 21"
Width between arms	49.5 cm / 20"
Weight of product	7.8 kgs / 17 lbs
Maximum user weight	133 kgs / 293 lbs

Optional Model no. 3521A
Adjustable Arms
Height: 57-72 cm / 22-28"

Optional Model no. 3521D
Detachable Arms

Optional Model no. 3521AD
Adjustable Height, Detachable Arms
Height: 54-69 cm / 21-27"



Stacking Commode

- Traditional armchair appearance
- Powder coated steel
- Hygienic pot (removable) with lid
- Upholstered backrest and overseat - flame retardant

Dimensions

Seat height	54.5 cm / 21.5"
Seat width	41 cm / 16"
Seat depth	48 cm / 19"
Width	57 cm / 22"
Height	81 cm / 32"
Depth	56 cm / 22"
Weight of product	7 kgs / 15 lbs
Maximum user weight	133 kgs / 293 lbs

Optional Model no. PL0068
European Pot and Lid



Folding Commode

- Folds easily for storage and transportation
- Fold back armrests
- Epoxy coated steel
- Hygienic pot (removable) with lid
- Upholstered backrest and overseat - flame retardant

Dimensions

Seat height	50 cm / 20"
Seat width	46 cm / 18"
Seat depth	45 cm / 18"
Width	51 cm / 20"
Height	90 cm / 35.5"
Depth	58.5 cm / 23"
Width between arms	46 cm / 18"
Weight of product	7 kgs / 15 lbs
Maximum user weight	133 kgs / 293 lbs

Optional Model no. PL0068
European Pot and Lid

