

The Proctogram Chair

Defecography Procedure Chair

Designed to improve
patient safety and
stability during
defecography
procedures.

Developed in conjunction with
radiologists, radiographers and
manual handling advisors in
the NHS.



Improved patient safety
and stability



Easier and safer
access



Braked Front Castors



Multi-functional
positioning



Retracting privacy
screen for dignity



Developed with NHS
imaging teams



Proctogram Procedure
Packs available for 3001
& 3002 models



A Solution Developed From A Clinical Need

Imaging teams highlighted the challenges of existing equipment used during defecography procedures, often lacking stability, accessibility and patient support.

Roma Medical worked in partnership with radiologists, radiographers and manual handling advisors within the NHS to design a chair that improves patient safety, dignity and procedural efficiency.

THE CHALLENGE

The original equipment lacked arm rails or integrated steps making access difficult and increasing falls risk

A crudely attached copper plate was also required to improve attenuation.



Key Features & Benefits

- 1



Safe to access Hand rail with folding footstep
- 2



Mobile for storage, 2 rear castors & front braked castors
- 3



Remote controlled forward / reverse aperture seat adjustment for precise imaging
- 4



Copper plates can be mounted on either side of the seat to improve attenuation in smaller rooms
- 5



Liquid fillable seat to reduce magnification
- 6



Retracting privacy/modesty screen

Built for Clinical Environments

- ✓ Designed for radiology and imaging departments
- ✓ Easy to clean, durable materials
- ✓ Supports infection prevention protocols
- ✓ Reliable, robust and built to last

Technical Specifications

Model No. 3001

Width Front 620mm / 24"	Floor to Step 1 150mm / 6"
Width Back 495mm / 19"	Step 1 to Step 2 150mm / 6"
Length 1100mm / 43"	Seat Width 450mm / 18"
Height (to top of handle) 1150mm / 45"	Seat Depth 450mm / 18"
Height (to top of backrest) 1115mm / 44"	Transit Footrest to Seat Top 390mm / 15"
Height (to top of seat) 730mm / 29"	Width Between Handles 555mm / 22"



Proctogram Procedure Pack

Clinical waste proctogram procedure packs available to purchase as a consumable. Suitable for 3001 & 3002 models. Please use code **3001/Bag** to order.

“Using high-quality, bespoke equipment has elevated our DPG service from a makeshift examination to a truly patient-centred radiological procedure. Patients undergoing these sensitive pelvic floor investigations deserve to be treated with the same dignity, respect and clinical excellence afforded to any other imaging examination. Innovations such as the Roma Proctogram chair and bespoke odourless waste bags help us deliver exactly that and set the standard of care we should be aspiring to deliver.

-Dr Adam Aboalkaz
Consultant Gastrointestinal, Genitourinary and Pelvic floor Radiologist