



Product Catalogue

Core Product Range

2026



MOVEMENT, MADE POSSIBLE.

At RMA, we believe movement is more than getting from A to B, it's about living life, uninterrupted.

For over 50 years, we've partnered with healthcare professionals, customers and communities to design mobility and care solutions that build confidence, restore independence and unlock potential.

Proudly manufacturing in the UK, everything we do is driven by insight, engineering for reliability and built to make a meaningful difference.

MORE THAN MOBILITY

Mobility • Clinical • Sport

CONTENTS

- 4 Boot Scooters
- 8 Mid Size Scooters
- 12 Road Scooters
- 16 Powered Wheelchairs
- 22 Manual Wheelchairs
- 28 Technical Specifications



BOOT SCOOTERS

COMPACT MOBILITY SCOOTERS

Our Boot Scooter mobility range incorporates the latest innovative designs for easy transportation to shops and leisure areas and are ideal for indoor use in confined spaces.

These Scooters can be easily disassembled without tools and will fit into the boot of most vehicles for easy transport.

Browse our selection to find the perfect scooter to meet your needs.

MORE THAN MOBILITY

FREEDOM

Mobility should fit around your life. Our products are selected to help make everyday journeys, days out and time with others more accessible.



Virgo
S110


M2089

TS30+


This scooter incorporates the latest innovative folding design and can be used both indoors and out, providing a practical mobility solution that is lightweight, compact and easy to transport. Designed to maximise your independence, it offers comfort, convenience and reliable everyday performance wherever your journey takes you.




Orange



Quick Folding Design




Comfortable Padded Seat




Lightweight Aluminium Frame



Airline-Friendly Lithium Battery



Easy to Transport



Solid Puncture Proof Tyres

SOLAX

Powered by a 25.2V lithium battery and 150W motor, the TS30+ delivers a top speed of 4mph and a range of up to 9 miles. Its quick-fold frame reduces to just 90 × 55 × 28cm, while solid puncture-proof tyres require minimal maintenance and a comfortable padded seat provides support throughout everyday journeys.

USER WEIGHT 125kg (19.7st)

4mph (6kmph)

9mi (14km)*

S110

Virgo

This scooter incorporates the latest innovative designs and can be used in confined spaces both indoors and out, aiming to provide you with maximum comfort whilst also being as safe and efficient as possible. This model comes with an interchangeable colour panel set including, Red, Grey and Blue panels.


- Flip Up Armrests
- Fold down Backrest
- Swivel Seat
- Rear Carry Handle
- Anti-Tipper Wheels
- Detachable Battery Pack
- Puncture Proof, Wheels
- Detachable Basket

USER WEIGHT

136kg (21.5st)

4mph (6kmph)


8mi (13km)*



Red



Grey



Blue

Interchangeable Colour Panels

**RMA**
MOBILITY

GK4

Alora

Say hello to the future of convenience with the SHOPRIDER® Alora! This cutting-edge mobility scooter is not just a mode of transport; it's a lifestyle upgrade. Available in a stunning Metallic Teal finish, the Alora is set to turn heads wherever you go.

- Horn
- Adjustable Armrest Width
- Fold down Backrest
- Adjustable Seat Height
- Rear Carry Handle
- Anti-Tipper Wheels
- Detachable Battery Pack
- Puncture Proof Wheels
- Detachable Basket

USER WEIGHT

114kg (18st)

4mph (6kmph)

8mi (13km)*



Teal

**SHOPRIDER®**

MID-SIZE SCOOTERS

STANDARD MOBILITY SCOOTERS

Our Mid-Size Mobility Scooters are designed for more substantial use than the mini boot scooters. They are still transportable and can travel over rougher terrain with an increased travel range and higher load.

Explore the mid-size scooter range to find the best scooter for your lifestyle.

MORE THAN MOBILITY

COMFORT

Every journey should feel comfortable and reassuring. That's why we consider the everyday details that help people travel with greater ease.




Valencia
GK10

GK9 Paris

The SHOPRIDER® Paris is a lightweight mid-size mobility scooter with a height adjustable padded seat & armrests that adjust and fold for easier entry. The tiller can also be adjusted to provide a comfortable amount of legroom.

- Fold down Backrest
- Flip Up Adjustable Armrests
- Adjustable Seat Height
- Adjustable Tiller Angle
- Detachable Basket
- Puncture Proof Tyres
- Non-Marking Wheels
- Detachable Battery Pack
- Anti-Tipper Wheels
- Rear Carry Handle

USER WEIGHT	136kg (21.5st)
	4mph (6kmph)
	16mi (24km)*



Charcoal



SHOPRIDER®

S888NR Sovereign

The Sovereign 4 is equipped with a high-performance electric motor that provides smooth and efficient navigation over steep hills and rough terrain with ease.

The scooter features a maximum speed of 4mph and up to 20 miles on a single charge, making it an excellent option for both short and long-distance travel.

- Front Shopping Basket
- Height & Depth Adjustable Seat
- Flip Up Armrests
- Fold Down Backrest
- Anti-Tipper Wheels
- Five Position Tiller Angle
- Chrome Bumper
- Headlights

USER WEIGHT	136kg (21.5st)
	4mph (6kmph)
	16mi (24km)*



Metallic Red



SHOPRIDER®

GK10 Valencia


The SHOPRIDER® Valencia Mobility Scooter has a range of up to 20 miles between charges and dismantles easily without tools for storage in a car boot. It is available in Blue or Red.




Metallic Blue



Metallic Red



Comfort can be tailored to the individual with a height-adjustable, four-position swivel seat, width-adjustable armrests and an adjustable-angle tiller. A headlight, rear reflectors and puncture-proof wheels provide added reassurance, while its 136kg maximum user weight makes the Valencia a comfortable and dependable choice for everyday mobility.



SHOPRIDER®

USER WEIGHT	136kg (21.5st)		4mph (6kmph)		20mi (32km)*
--------------------	-----------------------	---	---------------------	---	---------------------

ROAD SCOOTERS

8MPH ROAD MOBILITY SCOOTERS

Our 8mph Road Scooters are designed for heavy daily use and are capable of carrying larger weights over longer distances, generally alleviating the need for them to be transported.

Find the perfect scooter for your needs in our Road Scooter range.

**All class 3 scooters must be registered with the DVLA before use on the road*

MORE THAN MOBILITY

CONFIDENCE

Whether travelling locally or venturing further, dependable mobility can provide the confidence to enjoy greater freedom every day.



ES30x
R4SE

S9/R **Allrounder**


Introducing the Allrounder Road Scooter:
Versatile & Reliable. Features flip-up armrests, adjustable Delta tiller, and full suspension for comfort. With pneumatic tires, a shopping basket, and USB charging ports for convenience. Complete with LED lights, adjustable mirrors, and a nudge bumper for safety.

- Flip Up armrests
- Adjustable Delta Tiller (with storage)
- Full Suspension
- Black, Pneumatic, Non-Marking Tyres
- Shopping basket with carry handle
- USB 3.0 & USB C charging ports
- High Back, Contoured Seat
- Adjustable Rear View Mirrors
- LED Lights
- Emergency Brake

USER WEIGHT 159kg (25st)

8mph (12kmph)

20mi (32km)*



Red



SHOPRIDER®



R4SE **ES30x**

This scooter incorporates the latest luxury design and advanced technology to provide exceptional comfort and performance for everyday road use. Designed for longer journeys, it offers a smooth, stable and reliable ride, helping you maintain your independence with confidence wherever you travel.

- Heated Seat
- Heated Handlebars
- GPS Navigation
- Bluetooth Audio & Radio
- Reversing Camera
- Automatic Indicators & Headlights
- Pneumatic Tyres
- Up to 30 Mile Range

USER WEIGHT 160kg (26.7st)

8mph (12kmph)

30mi (45km)*



Black



RUIDI



P889SLB **Cadiz**


Introducing the SHOPRIDER® Cadiz:
a versatile mobility scooter with full suspension, adjustable seating, and width-adjustable arms for personalised comfort. Complete with a full lights set, high/low-speed switch, and front/rear suspension, it ensures a smooth and safe ride every time.

- High/Low Speed Switch
- Width Adjustable Arms
- Forward/Back/Height Adjustable Seat
- Full Suspension Front & Rear
- Large Pneumatic Tyres
- Full Lights Set
- Adjustable Angle Tiller with Delta Handle


USER WEIGHT 160kg (25st)

4 & 8mph (6 & 12kmph)

20mi (32km)*



Red



Black



SHOPRIDER®



889XLSBN/CAN **Traveso**


Introducing the SHOPRIDER® Traveso.
With headlights, indicators, and rear lights ensuring visibility, it offers safety wherever you go. Experience comfort with an adjustable captain seat inside a luxurious rust-proof cab.

- Headlights & Indicators
- High/Low Speed Switch
- Delta Tiller & Horn
- Windscreen Wiper & Washer
- Luxurious Rust Proof Cab
- Adjustable Captain seat with Headrest
- Window Vents
- Rear Lights
- Non-marking, Pneumatic Tyres
- Full Suspension
- Two Lockable Doors

USER WEIGHT 160kg (25st)

4 & 8mph (6 & 12kmph)


Up to 20mi (32km)*



Silver



SHOPRIDER®



POWER CHAIRS

POWER & ELECTRIC WHEELCHAIRS

Designed for both indoor and outdoor use, these chairs are built for users requiring a joystick control as opposed to a tiller control. Roma's Powerchairs are transportable and manoeuvrable in confined spaces.

Take a look through our catalogue of power chairs to best suit your everyday needs.

MORE THAN MOBILITY

INDEPENDENCE

The right mobility solution can make everyday activities more manageable, helping people maintain control, choice and independence.



Reno Elite
P325

P5519 Athens

Designed for nimble everyday mobility, the Athens is a lightweight powerchair that feels equally at home in compact indoor spaces and outdoors. Its aluminium frame, tight turning radius and responsive twin-motor drive make it particularly easy to manoeuvre, while the removable battery and portable design simplify transport and storage.

- Lightweight Aluminium Frame
- Compact 19.2kg Design
- Under 900mm Turning Radius
- Twin 250W Motors
- Removable Lithium Battery
- Up to 9.5 Mile Range
- Easy to Lift & Store
- Joystick Control

USER WEIGHT 136kg (21.5st)

4mph (6kmph)

9.5mi (15km)*

Black



RMA
MOBILITY

P600 EP30+

Created for users who want greater range without sacrificing portability, the EP30+ combines a lightweight folding frame with dual lithium batteries and dependable twin-motor performance. Offering up to 25 miles per charge and puncture-proof tyres, it is a practical choice for longer days out, everyday travel and straightforward storage.

- Aluminium & Magnesium Alloy Frame
- Quick Folding Design
- Dual Lithium Battery System
- Twin 250W Motors
- Solid Puncture-Proof Tyres
- Up to 25 Mile Range
- Built for Longer Journeys
- Joystick Control

USER WEIGHT 135kg (21.2st)

4mph (6kmph)

25mi (40km)*

Deep Green



RMA
MOBILITY

P5538 Rhodes

Built for convenient travel and everyday independence, the Rhodes combines a lightweight magnesium frame with a quick, compact fold. Strong yet easy to handle, it offers a comfortable and dependable drive indoors and out, with a 13-mile range that makes it well suited to daily journeys, car transport and storage.

- Lightweight Magnesium Frame
- Quick Compact Fold
- 19.9kg Including Battery
- Removable Lithium Battery
- 150W Twin Motors
- Up to 13 Mile Range
- Ideal for Car Transport
- Joystick Control

USER WEIGHT 120kg (19st)

4mph (6kmph)

13mi (21km)*

Black



RMA
MOBILITY

P325 Reno Elite

Introducing the Reno Elite Power Wheelchair. With a fold-down backrest, lap belt, and removable battery pack, it prioritises safety and convenience. Enjoy comfort with the adjustable footplate and padded seat, while the joystick controller offers easy operation.

- Fold Down Backrest
- Lap Belt
- Removable Battery Pack
- Puncture Proof, Non-Marking Wheels
- Flip Up & Height Adjustable Footplate
- Storage Basket
- Padded Seat
- Joystick Controller
- Width Adjustable, Flip Up Armrests

USER WEIGHT 136kg (21.5st)

4mph (6kmph)

12mi (20km)*

Flame Red

Flame Blue



RMA
MOBILITY

P325/CS

Reno Elite Captain Seat

Introducing the Reno Elite Captain Seat. Elevate your comfort and style with the luxurious captain seat. Available in Flame Red or Flame Blue, this chair is designed for ultimate relaxation pairing plush comfort with convenience. Make it your own with our many customisable features.

- Fold Down Backrest
- Lap Belt
- Removable Battery Pack
- Puncture Proof, Non-Marking Wheels
- Flip Up & Height Adjustable Footplate
- Storage Basket
- Padded Seat
- Joystick Controller
- Width Adjustable, Flip Up Armrests

USER WEIGHT 136kg (21.5st)

4mph (6kmph)

12mi (20km)*



Flame Red



Flame Blue



RMA
MOBILITY



P400

EP50+

This powerchair incorporates the latest innovative design to provide reliable indoor and outdoor mobility, offering a smooth, stable and comfortable driving experience. Designed to maximise your independence, it combines adjustable seating, powerful rear-wheel drive and everyday practicality to keep you moving with confidence.

- Rear-Wheel Drive
- Adjustable Seat & Backrest
- Dynamic Joystick Controller
- Electromagnetic Brakes
- Pneumatic Tyres
- Up to 28 Mile Range
- Up to 6 mph Speed
- 136kg Maximum User Weight

USER WEIGHT 136kg (21.5st)

6mph (10kmph)

28mi (49km)*



Orange



RMA
MOBILITY

P325/RS

Reno Seat Riser


Introducing the Reno Seat Riser Power Wheelchair. Revolutionising mobility with dynamic height adjustment. Sharing key features with the Reno Elite and Reno Elite Captain Seat, this model offers unparalleled versatility by adjusting its seat height on demand. Experience enhanced accessibility and convenience with customisable elevation.

- Fold Down Backrest
- Removable Battery Pack
- Puncture Proof, Non-Marking Wheels
- Flip Up & Height Adjustable Footplate
- Padded Seat
- Joystick Controller
- Width Adjustable, Flip Up Armrests

USER WEIGHT 136kg (21.5st)

4mph (6kmph)

12mi (20km)*



Flame Red



Flame Blue



RMA
MOBILITY

P400 /RS

EP50+ Riser

This powerchair incorporates an innovative 300mm seat riser, providing greater reach and flexibility throughout the day. Designed to maximise your independence, it combines adjustable seating, powerful rear-wheel drive and everyday practicality for smooth, stable and comfortable indoor and outdoor mobility.

- 300mm Seat Riser
- Rear-Wheel Drive
- Adjustable Seat & Backrest
- Dynamic Joystick Controller
- Electromagnetic Brakes
- Pneumatic Tyres

USER WEIGHT 136kg (21.5st)

6mph (10kmph)

28mi (49km)*



Orange



RMA
MOBILITY

MANUAL WHEELCHAIRS TRANSIT & SELF PROPELLED

Designed for both indoor and outdoor use, Roma's manual wheelchairs are built to provide comfort, practicality and dependable everyday mobility. With a range of self-propelled and transit models available, there's a wheelchair to suit a variety of users and lifestyles.

Take a look through our catalogue of manual wheelchairs to find the model that best suits your everyday needs.

MORE THAN MOBILITY

SUPPORT

With over 50 years of mobility experience, RMA is committed to providing dependable products and knowledgeable support when it matters.



Wisp
1650

1650 / 1640

Wisp

This wheelchair incorporates the latest lightweight design to provide comfortable and reliable mobility for everyday use. Designed to maximise your independence, it combines a lightweight folding frame with practical features for easy transportation, storage and everyday convenience wherever you go.

- Lightweight Aluminium Frame
- Quick Release Rear Wheels
- Half Folding Backrest
- Flip-Up Removable Armrests
- Attendant Cable Brakes
- Removable Leg-rests


RMA
MOBILITY


	SP	CT
USER WEIGHT	18.8st (120kg)	19st (120kg)
KG	12.5kg	11.5kg



Orange



Grey



Green



Blue

1610 / 1630

Avant

This wheelchair combines durable steel construction with practical everyday features to provide comfortable and reliable mobility. Designed to maximise your independence, it offers a smooth self-propelling experience whilst remaining easy to transport and store thanks to its folding design.

- Strong Steel Frame
- Fold Down Backrest
- Detachable Desk Armrests
- Swing Away Footrests
- Crash Tested


RMA
MOBILITY

	SP	CT
USER WEIGHT	17.9st (114kg)	17.9st (114kg)
KG	17.5kg	17.5kg



Silver

1661 / 1660

Lightweight Wheelchair

This wheelchair combines a lightweight aluminium frame with practical everyday features to provide comfortable and reliable mobility. Designed to maximise your independence, it is easy to transport, store and use, making it an ideal choice for everyday travel and daily use.

- Lightweight Aluminium Frame
- Quick Release Rear Wheels
- Flip-Out & Removable Leg Rests
- Fold Down Handles
- 120kg (18.9st) Maximum User Weight


RMA
MOBILITY

	SP	CT
USER WEIGHT	18.8st (120kg)	18.9st (120kg)
KG	10.2kg	9.2kg



Light Blue



Silver

P001

Power Pack

The Roma Power Pack is engineered to aid in the pushing of the wheelchair, delivering additional assistance on inclines and overall support for the carer. Features an easy-to-operate hand control lever and security switch key for safety. With twin drive wheels and a removable battery pack, it offers enhanced mobility and even easily dismantles for transportation and storage convenience.

- Easy to Operate Hand Control Lever
- Security Switch Key
- Twin Drive Wheels
- Removable Battery Pack Remote Charging
- Dismantles Easily for Transportation & Storage

USER WEIGHT	136kg (21.5st)
	3.4mph (5.5kmph)
	10mi (16km)*


RMA
MOBILITY

SHOPRIDER® SCOOTERS

A TRUSTED NAME IN MOBILITY

At Roma Medical we pride ourselves on offering the pinnacle of mobility solutions. As an esteemed retailer of SHOPRIDER® products, we stand by their unparalleled quality and reliability.

Each SHOPRIDER® scooter is meticulously crafted to adhere to stringent ISO9001 standards, ensuring exceptional durability and performance. Our commitment to excellence extends beyond mere functionality; it encompasses innovative design features that enhance user experience and promote independence.

With years of experience in the market, our customers can rest assured knowing they are investing in a product that is built to last.

SHOPRIDER®



AllRounder
S9/R

Model Photo			
Model Name	TS30+	Alora	Virgo
Model No.	M2089	GK4	S110
Class	2	2	2
Battery Capacity	25.2V 10Ah Lithium	2 x 12V 12Ah	2 x 12V 12Ah
Max. User Weight	125kg (19.7st)	114kg (18st)	136kg (21.5st)
Range per Charge*	9mi (14.4km)	8mi (13km)	8mi (13km)
Turning Radius	1200mm (47")	1620mm (63")	1150mm (43.3")
Max. Speed	4mph (6kmph)	4mph (6kmph)	4mph (6kmph)
Safe Working Angle	6° (10%)	6° (10%)	6° (10%)
Height	960mm (37") Folded: 280mm (11")	850mm (33")	1100mm (43.3")
Width	550mm (21.6") Folded: 550mm (21.6")	510mm (20")	560mm (22")
Length	990mm (39") Folded: 900mm (35.4")	1020mm (40")	1050mm (41.3")
Clearance	30mm (1.2")	50mm (2")	65mm (2.6")
Total Weight	19.6kg (43lb)	40kg (90lb)	43.5kg (96lb)
Colours	Orange	Teal	Red/Grey/Blue Interchangeable Colour Panels

*Passenger weight, way of driving, terrain and ground conditions, temperatures, weather conditions and battery conditions can affect maximum range.

Model Photo			
Model Name	Paris	Sovereign	Valencia
Model No.	GK9	S888NR	GK10
Class	2	2	2
Battery Capacity	2 x 21V 12Ah	2 x 12V 34/36Ah	2 x 12V 34/36Ah
Max. User Weight	136kg (21.5st)	136kg (21.5st)	136kg (21.5st)
Range per Charge*	16mi (25km)	20mi (32km)	20mi (32km)
Turning Radius	1150mm (45")	1150mm (45")	1330mm (52")
Max. Speed	4mph (6kmph)	4mph (6kmph)	4mph (6kmph)
Safe Working Angle	8° (14%)	8° (14%)	8° (14%)
Height	880mm (35")	940mm (37")	960mm (38")
Width	560mm (22")	580mm (23")	585mm (23")
Length	1050mm (41")	1250mm (49")	1170 mm (46")
Clearance	75mm (3")	75mm (3")	90mm (3.5")
Total Weight	56kg (120lbs)	82kg (180lbs)	75kg (165.5lbs)
Colours	Charcoal/Burgundy	Metallic Red	Metallic Red/Blue

Model Photo				
Model Name	Allrounder	Cadiz	ES30x	Traveso
Model No.	S9/R	P889SLB	R4SE	S889XLSBN/CAN
Class	3	3	3	3
Battery Capacity	2 x 12V 34/36Ah	2 x 12V 34/36Ah	2 x 24V 75Ah	2 x 12V 36/75Ah
Max. User Weight	159kg (25st)	160kg (25st)	160kg (25st)	160kg (25st)
Range per Charge*	20mi (32km)	20mi (32km)	30 mi (48km)	20mi (32km)
Turning Radius	1500mm (59")	1500mm (59")	2100mm (82")	1780mm (70")
Max. Speed	8mph (12kmph)	4/8mph (6/12kmph)	8mph (12kmph)	4/8mph (6/12kmph)
Safe Working Angle	9° (15%)	8° (14%)	12° (21%)	8° (14%)
Height	1300mm (51")	1120mm (44")	1300mm (51")	1620mm (63.7")
Width	630mm (24")	640mm (25")	720mm (28")	730mm (28.7")
Length	1365mm (53")	1300mm (51")	1420mm (56")	1600 mm (63")
Clearance	80mm (3")	150mm (5.9")	120mm (4.7")	70mm (2.7")
Total Weight	107kg (235lbs)	100kg (220lbs)	157kg (346lbs)	150kg (330lbs)
Colours	Red	Red/Black	Black	Silver

*Passenger weight, way of driving, terrain and ground conditions, temperatures, weather conditions and battery conditions can affect maximum range.
**Safe working angle per user weight; 0° (0%) - 225kg (35 st 6lbs) / 6° (10%) - 150kg (23st 8lbs) / 9° (15%) - 125kg (19st 9lbs) / 12° (20%) - 100kg (15st 10lbs)
***Class 3 scooters must be registered with the DVLA

Model Photo							
Model Name	Wisp		Lightweight Wheelchair		Avant		Power Pack
Model No.	1650	1640	1661	1660	1610	1630	P001
Type	SP	CT	SP	CT	SP	CT	Power Assistance
Overall Width	670mm	670mm	690mm	630mm	660mm	660mm	200mm
Overall Length	990mm	915mm	1000mm	1000mm	1070mm	1070mm	245mm
Overall Height	950mm	950mm	920mm	920mm	910mm	910mm	300mm
Product Weight	12.5kg	11.5kg	10.2kg	9.2kg	17.5kg	17.5kg	14kg
Max. User Weight	120kg	120kg	120kg	120kg	114kg	114kg	136kg
Colours	Blue Green Grey Orange	Blue Green Grey Orange	Silver Blue	Silver Blue	Silver	Silver	Black

Model Photo				
Model Name	Athens	Rhodes	EP30+	Reno Elite
Model No.	P5519	P5538	P600	P325
Class	2	2	2	2
Battery Capacity	2 x 24V 10Ah	2 x 24V 12Ah	2 x 24V 10.6Ah	2 x 12V 22Ah
Max. User Weight	136kg (21.5st)	120kg (18.9kg)	135kg (21.2st)	136kg (21.5st)
Range per Charge*	9.5mi (15km)	13mi (21km)	25mi (40km)	12mi (20km)
Turning Radius	<900mm (35.5")	<1200mm (47.2")	<1200mm (47.2")	550mm (21.6")
Max. Speed	4mph (6kmph)	4mp (6kmph)	4mph (6kmph)	4mph (6kmph)
Safe Working Angle	6° (10%)	6° (10%)	6° (10%)	6° (10%)
Height	480mm (18.8")	1030mm (40.5")	1060mm (41.7")	950mm (37")
Width	440mm (17.3")	605mm (23.8")	600mm (23.6")	600mm (24")
Length	400mm (15.7")	940mm (37")	900mm (35.5")	880 mm (35")
Clearance	40mm (1.5")	40mm (1.5")	40mm (1.5")	50mm (2")
Total Weight	17.7kg (39lbs)	19.9kg (43.8lbs)	21.5kg (47lbs)	62kg (134lbs)
Colours	Black	Black	Green	Flame Red/Blue

*Passenger weight, way of driving, terrain and ground conditions, temperatures, weather conditions and battery conditions can affect maximum range.
**Safe working angle per user weight; 0° (0%) - 225kg (35 st 6lbs) / 6° (10%) - 150kg (23st 8lbs) / 9° (15%) - 125kg (19st 9lbs) / 12° (20%) - 100kg (15st 10lbs)
***Class 3 scooters must be registered with the DVLA

Model Photo				
Model Name	Reno Elite CS	Reno Seat Riser	EP50+	EP50+ Riser
Model No.	P325/CS	P325/RS	P400	P400 Riser
Class	2	2	3	3
Battery Capacity	2 x 12V 22Ah	2 x 12V 22Ah	2 x 24V 52Ah	2 x 24V 52Ah
Max. User Weight	136kg (21.5st)	136kg (21.5st)	136kg (21.5st)	136kg (21.5st)
Range per Charge*	12mi (20km)	12mi (20km)	28mi (45km)	28mi (45km)
Turning Radius	550mm (21.6")	550mm (21.6")	900mm (35.5")	900mm (35.5")
Max. Speed	4mph (6kmph)	4mph (6kmph)	6mph (10kmph)	6mph (10kmph)
Safe Working Angle	6° (10%)	6° (10%)	10° (16.67%)	10° (16.67%)
Height	1000mm (41")	1130-1410mm (44-55")	1215mm (48")	1350mm - 1650mm (53" - 65")
Width	600mm (23.6")	600mm (24")	580mm (23")	580mm (23")
Length	880mm (34.6")	880mm (35")	1080mm (42.5")	1080mm (42.5")
Clearance	50mm (2")	50mm (2")	50mm (2")	50mm (2")
Total Weight	66.5kg (146.6lbs)	79.5kg (175lbs)	106.5kg (235lbs)	136kg (299lbs)
Colours	Flame Red/Blue	Flame Red/Blue	Orange	Orange

MORE THAN MOBILITY

MOBILITY

CLINICAL

SPORT

Why should I buy an RMA product?



EXPERIENCE

RMA has decades of experience in mobility. Since 1970, we have been designing, building and supplying a variety of products that work in the real world.



CHOICE

As this catalogue shows, RMA is able to offer real choice and real options, with a broad range covering mobility, daily living and specialist products.



SUPPORT

We don't disappear once the order has been delivered. We are in it for the long haul: developing long-term products and building long-term relationships. Our team is here to support you and your customers.



VALUE

All of us want good value for money. RMA products are designed and built to offer the right balance of quality, performance and price.



INNOVATION

As engineers and designers, we are continually developing our products to make them better, easier to use and more relevant. If you have an idea for a better way to do something, please share it with us.



PEOPLE

Many of our colleagues have been working for RMA for over 20 years. We know that behind every product is a team that keeps the needs of the customer at the centre of everything it does.



WHAT TO CONSIDER

Thinking through the complete journey can help narrow down the most suitable option.

1. HOW WILL YOU USE IT?




Consider where you plan to travel, how often it will be used and the surfaces you are likely to encounter.

2. COMFORT & CONTROL




Check that the seating, legroom, support & controls feel comfortable and can be operated without strain.

3. STORAGE & CHARGING



Measure doorways, gates and storage areas. Check that a suitable, dry charging location is available.

4. TRANSPORT & LIFTING




Check the folded dimensions, total weight and the heaviest individual component, and whether it will fit the intended vehicle.

5. RANGE & CAPACITY



Choose a range and maximum user weight suited to typical journeys. Hills, temperature & additional loads may affect performance.

6. TRY IT



Whenever possible, try the product. Ask about demonstrations, warranty, servicing, replacement parts and dealer support.




Think through the complete journey, from leaving home and loading the product to charging and storing afterwards.

YOUR MOBILITY DEALER CAN HELP

This guide is a starting point. Your local mobility dealer can explain the differences, check suitability and arrange a demonstration where available.

Requirements may vary depending on the product, intended use and transport provider. Your mobility dealer can provide further guidance.





E: sales@romamedical.co.uk
W: www.romamedical.co.uk
T: 01656674488

York Road, Bridgend Industrial Estate,
Bridgend, CF31 3TB

SHOPRIDER®